

ZHERMACK S.P.A.

Safety Data Sheet according to the Dir. 2001/58/EEC

GINGIFAST ELASTIC

Revision date: 09/06/03

Print Date: 03/08/2007

1. IDENTIFICATION OF THE PREPARATION AND OF THE COMPANY

1.1 Identification of the preparation:

Trade name: GINGIFAST ELASTIC

Trade code: C203225 (single pack)/C401500 (blister)

1.2 Utilization of the substance or of the preparation:

Gingival of polivinilsiloxane for gingival masks reproduction.

1.3 Identification of the company:

Supplier: ZHERMACK S.P.A. VIA Bovazecchino, 100 -45021- Badia Polesine (RO)
ITALY

Emergency telephone number of the company: 0039 0425 597611

2. COMPOSITION/INFORMATION ON INGREDIENTS

Hazardous components within the meaning of EEC directive 67/548 and corresponding classification: None

3. HAZARDS IDENTIFICATION

PREPARATION NOT CLASSIFIED DANGEROUS

No specific hazards are encountered under normal product use.

4. FIRST AID MEASURES

4.1 Contact with skin: Wash with plenty of water and soap.

4.2 Contact with eyes: Wash immediately with water for at least 10 minutes.

4.3 Swallowing: Induce vomiting. SEEK A MEDICAL EXAMINATION IMMEDIATELY and present the safety-data sheet. A suspension of activated charcoal in water, or liquid paraffin may be administered.

4.4 Inhalation: Ventilate the premises.

5. FIRE-FIGHTING MEASURES

5.1 Recommended extinguishers: water, CO₂, foam, chemical powders, according to the materials involved in the fire.

5.2 Extinguishers not to be used: none in particular.

5.3 Risks arising from combustion: avoid inhaling the fumes.

5.4 Protective equipment: none in particular.

6. ACCIDENTAL RELEASE MEASURES

6.1 Measures for personal safety: Use gloves and protective clothing.

6.2 Environmental measures: limit leakages with earth or sand. If the product has escaped into a water course, into the drainage system, or has contaminated the ground or vegetation, notify the competent authorities.

| REVISIONE | DATA | DESCRIZIONE | EMISSIONE | APPROVAZIONE |
|-----------|----------|--------------------------------|-----------|--------------|
| 1 | 01/03/99 | Prima emissione | | |
| 2 | 09/06/03 | Aggiornamento Dir. 2001/58/CEE | R&S | RSGQA |

ZHERMACK S.P.A.

Safety Data Sheet according to the Dir. 2001/58/EEC

GINGIFAST ELASTIC

Revision date: 09/06/03

Print Date: 03/08/2007

6.3 **Cleaning methods:** recover the product for re-use if possible, or for elimination. After the product has been recovered, rinse the area and materials involved with water.

7. HANDLING AND STORAGE

7.1 **Handling precautions:** See, too, paragraph 8 below. Do not eat or drink while working.

7.2 **Incompatible materials:** none in particular. See, too, paragraph 10 below.

7.3 **Storage conditions:** keep the closed container a temperature varying from 5°C to 27°C (41–80°F).

7.4 **Technical measures:** satisfactory room ventilation.

8. EXPOSURE CONTROLS / PERSONAL PROTECTION

8.1 **Precautionary measures:** give adequate ventilation to the premises where the product is stored and/or handled.

8.2 **Respiratory protection:** not needed for normal use

8.3 **Protection for hands:** not needed for normal use

8.4 **Eye protection:** not needed for normal use

8.5 **Protection for skin:** no special precaution must be adopted for normal use.

8.6 **Exposure limit(s) (ACGIH):** none

9. PHYSICAL AND CHEMICAL PROPERTIES

9.1 **Appearance and colour:** pink fluid vascularized.

9.2 **Odour:** odourless

9.3 **Relative density:** 1.125 g/cc

9.4 **Solubility in water:** insoluble

10. STABILITY AND REACTIVITY

10.1 **Conditions to avoid:** stable under normal conditions.

10.2 **Substances to avoid:** none in particular.

10.3 **Hazardous decomposition products:** none

11. TOXICOLOGICAL INFORMATION

The product does not contain toxicologically relevant substances. See, too, paragraph 2.

12. ECOLOGICAL INFORMATION

Adopt sound working practices, so that the product is not released into the environment.

13. DISPOSAL CONSIDERATIONS

Recover if possible. In so doing, comply with the local and national regulations currently in force.

14. TRANSPORT INFORMATION

Not classified as dangerous in the meaning of transport regulations.

15. REGULATORY INFORMATION

Page 2 of 2

| REVISIONE | DATA | DESCRIZIONE | EMISSIONE | APPROVAZIONE |
|-----------|----------|--------------------------------|-----------|--------------|
| 1 | 01/03/99 | Prima emissione | | |
| 2 | 09/06/03 | Aggiornamento Dir. 2001/58/CEE | R&S | RSGQA |

ZHERMACK S.P.A.

Safety Data Sheet according to the Dir. 2001/58/EEC

GINGIFAST ELASTIC

Revision date: 09/06/03

Print Date: 03/08/2007

The preparation should not be considered as dangerous according to dir. 88/379/EEC.

16. OTHER INFORMATION

The information contained herein is based on our state of knowledge at the above-specified date. It refers solely to the product indicated and constitutes no guarantee of particular quality. It is the duty of the user to ensure that this information is appropriate and complete with respect to the specific use intended. This MSDS cancels and replaces any preceding release.

Page 3 of 2

| REVISIONE | DATA | DESCRIZIONE | EMISSIONE | APPROVAZIONE |
|-----------|----------|--------------------------------|-----------|--------------|
| 1 | 01/03/99 | Prima emissione | | |
| 2 | 09/06/03 | Aggiornamento Dir. 2001/58/CEE | R&S | RSGQA |