



MGUARD

MGUARD DISPOSABLE ISOLATION GOWN



MGuard Gown Blue Level 1 Protection

Product Code: 70110

Carton of 50

Light weight gowns are designed for procedures where there is low exposure to aerosol.

- Level 1 Barrier Protection
- Light Weight
- Cool and Comfortable
- Splash resistant
- Non-woven
- Stretch cotton cuffs
- Neck and waist ties
- Universal size



MGuard Gown Medium Level 2 Protection

Product Code: 70120

Carton of 20

Medium weight gown where there is a risk of large splashes of blood or body substances.

- Level 2 Barrier Protection
- Medium Weight
- Stretch cotton cuffs
- Non-woven
- Neck and waste ties
- Universal size

ISOLATION CLOTHING FABRIC



Non-woven fabric with film

One way ventilation

High-efficiency anti bacteria

Thickened non-woven material through hot fusion and machine sewing, comfortable, breathable and anti-fouling. Oil and water repellent film prevents cross-infection and keeps personal hygiene clean.



Neck Tie



Ultrasonic Sewing



Knitted Cuff